

Technical Note

Revision history for single grain uniformity correction software: CorrSGbin

Version	Description	Date/ Responsible
1.0.0	First version of the program for correcting .bin files	2012/TLAP
1.1.0	Updated version that also allow correction of .binx files	11/2-14/TLAP
1.1.1	V.1.1.0 had a serious bug for correcting .binx files. This has been corrected and a check for availability of coordinate data for all grain positions (1-100) has been added	21/3-14/TLAP
1.2	This version also supports V7 bin files	22/7-15/TLAP
1.3	This version also supports V8 bin files	14/3-16/TLAP
1.4	Fixed "Range check error" resulting in clipped files.	2/5-18/KARB
1.5	Added some statistics for internal Risø use.	18/6-18/KARB